

Economical Dual Axis Power

A medium capacity, economical version of Harmar’s AL425 2-Axis lift for transporting power chairs and scooters. Double axis lift offers powered rotation and lifting, with adjustable arm height.

Features & Benefits

- Dual axis for powered rotation and lifting
- Adjustable arm height
- Weather-resistant hand control; wireless option available
- Passenger side standard; driver side available
- Offset post allows the mobility device to be safely unloaded on the side of the vehicle
- Docking device included based on application



Specifications	AL225 Axis II Lite Lift
Accommodates	Power chairs & scooters
Lifting Capacity	250 lbs. (113 kg)
Installed Weight	55 lbs. (25 kg)
Adjustable Boom Height	30" – 41"
Adjustable Boom Length	21.25" – 32.35"
Adjustable Boom Pitch	5°
Securement	Docking device included
Options	Battery pack, wireless remote
Warranty	3-year transferable

Key Features



Easy-to-use hand control with Soft-Touch buttons.



Convenient wireless key fobs available